

Name _____

Date _____

DNA Extraction Worksheet: Strawberry

Background: The long, thick fibers of DNA store the information for the functioning of the chemistry of life. DNA is present in every cell of plants and animals. The DNA found in strawberry cells can be extracted using common, everyday materials. We will use an extraction buffer containing salt, to break up protein chains that bind around the nucleic acids, and dish soap to dissolve the lipid (fat) part of the strawberry cell wall and nuclear membrane. This extraction buffer will help provide us access to the DNA inside the cells.


Pre-lab questions:

What do you think the DNA will look like?

Where is DNA found?

Materials:

- heavy duty ziploc bag
- 1 strawberry
- 10 ml DNA extraction buffer (soapy, salty water)
- cheesecloth
- funnel
- 50ml vial / test tube
- glass rod, inoculating loop, or popsicle stick
- 20 ml ethanol



Procedure:

1. Place one strawberry in a Ziploc bag.
2. Smash/grind up the strawberry using your fist and fingers for 2 minutes. *Careful not to break the bag!!*
3. Add the provided 10mL of extraction buffer (salt and soap solution) to the bag.
4. Knead/mush the strawberry in the bag again for 1 minute.
5. Assemble your filtration apparatus as shown to the right.
6. Pour the strawberry slurry into the filtration apparatus and let it drip directly into your test tube.
7. Slowly pour cold ethanol into the tube. OBSERVE _
8. Dip the loop or glass rod into the tube where the strawberry extract and ethanol layers come into contact with each other. OBSERVE _

Conclusions and Analysis

1. It is important that you understand the steps in the extraction procedure and why each step was necessary. Each step in the procedure aided in isolating the DNA from other cellular materials. Match the procedure with its function:

PROCEDURE

- A. Filter strawberry slurry through cheesecloth
- B. Mash strawberry with salty/soapy solution
- C. Initial smashing and grinding of strawberry
- D. Addition of ethanol to filtered extract

FUNCTION

- ___ To precipitate DNA from solution
- ___ Separate components of the cell
- ___ Break open the cells
- ___ Break up proteins and dissolve cell Membranes

2. What did the DNA look like? Relate what you know about the chemical structure of DNA to what you observed today.

3. Explain what happened in the final step when you added ethanol to your strawberry extract. (*Hint: DNA is soluble in water, but not in ethanol*)

4. A person cannot see single cotton thread 100 feet away, but if you wound thousands of threads together into a rope, it would be visible much further away. Is this statement analogous to our DNA extraction? Explain.

5. Why is it important for scientists to be able to remove DNA from an organism? List two reasons.

6. Is there DNA in your food? _____ how do you know?